


myLDR.com

Advanced Online Account Management

With **myLDR.com** you can

- Efficiently manage and maintain your dosimetry program online: check unreturned badges, browse dose history, export data into Excel® spreadsheets
- Give individual users read access to personal dose information by viewing the **Individual Dose Report**
- **Be GREEN** - go paperless while accessing all dose reports, invoices, E sign annual dose reports etc. online

A CLOSER LOOK...




myLDR Menu

Quick Add Button
for people & dosimeters



myLDR Landing Page

myLDR Dashboards




Receive critical alerts, e.g.,
ALARA & OVER Exposures

Track dosimeter shipments to
and from LANDAUER

Add, change, remove or view subaccounts,
participants and badges

View and create Dosimetry
Reports with Real Time Data

Access or download critical dose,
unreturned dosimeters or other data for
further analysis



GET STARTED WITH IDR TODAY...

■ Request an Individual Dose Report (IDR) login for your account.

This is a generic login that will be distributed to all participants.

Participants do **not** have the ability to change the password for the account login or access any functionality within myLDR.com.

■ Distribute the login to all badged employees

Use the login screen through myLDR.com to launch LDR

■ Individuals can enter the account & serial number listed on the back of the dosimeter.

■ Once logged in, you can view and export the most recent dose information and dose detail

They will only see information relating to their personal dose

End Session

Landauer Individual Dose Report

Please enter Account Number and Serial Number from your dosimeter and click the "Submit" Button.

Individual Dose Report (IDR) - System Demonstration

Individual Dose Report (IDR) - Glossary

Account Number *

Serial number *

* Required



To protect your privacy, no personal information is displayed

History Results

The doses are displayed in mSv

Dose Results

Total Records: 4

Doses as of 2012/10/31 10:47 CST

Dose Period	Total DDE	Total LDE	Total SDE	Beta	Total Neutron	Extremity
Q42012						M
2012						M
2011						
Lifetime	M	M	M			M

Visit myLDR.com to request your login today

Questions? Contact myLDRSupport@landauer.com or +44 (0) 1865 373008